


BROEN-LAB UniFlex™ CONNECTIONS

FOR TECHNICAL AND BURNING GAS

(Type F acc. to Doc. no.: 99G0020)



Benyt kun godkendte fittings. F.eks. VSH kompressionsfittings
Only use approved fittings. E.g. VSH compression fittings
Nur zugelassene Verschraubungen benutzen. Zum
Beispiel VSH Kompressionsverschraubungen

www.vsh-fittings.com

BROEN-LAB

www.BROEN-LAB.com CFR03022014 94G0153